




USB to RS-232 Adapter

Skye's DataHog dataloggers have a RS232 serial communications output. New PCs and laptops are often being supplied with only USB ports with no serial connections.

An alternative is to use a good quality, USB-Serial converter like the USB2-F-1001-A. Skye have tested this model with Windows XP, Vista, Windows 7 & 8. It is also suitable for use with Skye's HydroSense and SpectroSense meters.



Please note that some USB converters from other suppliers may not work with Skye dataloggers and meters. The USB converter isn't supplied with any software, so please go to the following website to download the relevant drivers.

<http://www.easysync-ltd.com/product/651/usb2-f-1001-a.html>

Also please click the following link for changing comm. ports when using these USB adaptors.

Version 1.2

Skye Instruments Ltd

21, Ddole Enterprise Park, Llandrindod Wells
Powys LD1 6DF, United Kingdom

TEL +44 (0)1597 824811

EMAIL skyeemail@skyeinstruments.com

WEB <http://www.skyeinstruments.com>